

# Scanning Vibrating Electrode Technique Time-Lapse Imaging (SVET-TLI): Studying Sacrificial Zinc Coatings in Salt Solution

**Rebecca Bolton**

Academic Supervisors: Prof. Geraint Williams, Prof. James Sullivan

Industrial Supervisor: Dr. Jörgen van de Langkruis

M2A has been made possible by the EU's Convergence European Social Fund through the Welsh Government



Swansea University  
Prifysgol Abertawe

**EPSRC**


Engineering and Physical Sciences  
Research Council




EPSRC Centre for Doctoral Training in Industrial Functional Coatings




Cronfa Gymdeithasol Ewrop  
European Social Fund




Scanning Vibrating Electrode Technique




Time Lapse Imaging




SVET-TLI



Results



Conclusions



Scanning Vibrating Electrode Technique




Time Lapse Imaging



SVET-TLI





Results



Conclusions

# Scanning Vibrating Electrode Technique


- In-situ measurement of localised corrosion activity
- Detects potential gradients produced by local currents generated by actively corroding surfaces



SVET probe assembly, adapted from Williams and McMurray, 2008 [1]

# Scanning Vibrating Electrode Technique

- Produces time-based, spatially-resolved current density maps



Red = Anodic = Metal Dissolution

Blue = Cathodic = Oxygen Reduction

Typical SVET map for Magzinc in 5% w/v NaCl after 24 hours

# Scanning Vibrating Electrode Technique

## Advantages


- Relatively non-perturbing
- Accurate qualitative and quantitative data that is both spatially and temporally resolved

## Disadvantages:

- Resolution restricted to probe diameter (100 $\mu$ m)
- No information on visual/physical behaviour



Scanning Vibrating Electrode Technique




Time Lapse Imaging



SVET-TLI




Results




Conclusions

- Images captured periodically to monitor activity or behaviour
- Regular intervals
- Automated software allows for long duration studies


# Time Lapse Imaging (in situ)




(a)



(b)



(c)



(d)

Photographic images of magnesium freely corroding in aerated 5% w/v NaCl (aq) taken at (a) 4, (b) 16, (c) 27, and (d) 38 min following immersion [1]

## Advantages:

- Potentially high resolution
- Generates progressive visual data
- Visual progression can also be used for rate data

## Disadvantages:


- No information other than visual




Scanning Vibrating Electrode Technique




Time Lapse Imaging



SVET-TLI





Results



Conclusions


- Combination of both SVET and time lapse photography
- 3D printed friction clamp
- Rest setting: probe assembly moves in the x direction to line up microscope
- Snapshot taken every minute






Scan Mode


Electrolyte



(a)



(b)



(c)




Scanning Vibrating Electrode Technique



Time Lapse Imaging



SVET-TLI




Results



Conclusions


- Zn coatings:
  - Hot Dip Galvanised (HDG)
  - Electrogalvanised (EG)
  - Physical Vapour Deposited (PVD0)
- Electrolytes:
  - 1% w/v NaCl
  - 5% w/v NaCl
- 6mm x 6mm surface scan area




HDG



EG




PVD0




# Results

## 1% NaCl

- Hot-Dip Galvanised Zinc







# Results

## 5% NaCl

- Hot-Dip Galvanised Zinc







# Results

## 1% NaCl

- Electrogalvanised Zinc







# Results

## 5% NaCl

- Electrogalvanised Zinc






# Results


## 1% NaCl

- Physical Vapour Deposited Zinc



00:00:00






# Results

## 5% NaCl

- Physical Vapour Deposited Zinc



00:00:00







Scanning Vibrating Electrode Technique




Time Lapse Imaging



SVET-TLI



Results





Conclusions

# Conclusions


	Hot Dip Galvanised	Electrogalvanised	Physical Vapour Deposited
1% NaCl	<ul style="list-style-type: none"><li>• Approx. 30 total anodes</li><li>• Rapid growth to a max size (600<math>\mu</math>m)</li><li>• Growth restricted by corrosion product</li></ul>	<ul style="list-style-type: none"><li>• 150+ anodes</li><li>• Rapid nucleation (200<math>\mu</math>m) and passivation followed by reactivation (600<math>\mu</math>m+)</li><li>• Entire surface corroded</li></ul>	<ul style="list-style-type: none"><li>• Approx. 100 anodes</li><li>• Rapid nucleation (400<math>\mu</math>m) followed by agglomeration</li><li>• Entire surface corroded</li></ul>
5% NaCl	<ul style="list-style-type: none"><li>• One large anodic front</li><li>• Corrosion product deposited in wake of the front</li><li>• Anodic activity restricted by corrosion product</li></ul>	<ul style="list-style-type: none"><li>• One large anodic front</li><li>• Corrosion product deposited either side and in wake of the front</li><li>• Anodic activity restricted by corrosion product</li></ul>	<ul style="list-style-type: none"><li>• Three large, expanding anodes</li><li>• Corrosion product deposited in wake of expansion</li><li>• Fronts meet and cease anodic activity</li></ul>

# Conclusions


- Photographic images display the visual effects of corrosion activity:




1 Hour



2 Hours



3 Hours




4 Hours

HDG Zinc in 5% NaCl


# Conclusions

- Photographic images resolve areas of net activity in SVET maps









Electrodeposited Zinc  
1% w/v NaCl  
At 1 hour

Electrodeposited Zinc  
1% w/v NaCl  
At 8 hours



# Conclusions

- Photographic images complement SVET data for comparative studies



PVD Zinc in 1% NaCl at 6hrs and 12hrs

PVD Zinc in 5% NaCl at 6hrs and 12hrs

- SVET map anodic and cathodic regions correlate with visual corrosion behaviour in photographs
- SVET-TLI is a useful technique to monitor both electrochemical and visual behaviour simultaneously and automatically
- Results show variation in corrosion mechanisms due to both the difference in coating production and electrolyte concentration

**TATA STEEL**

**M<sup>2</sup>A**  
MATERIALS AND MANUFACTURING ACADEMY

Thank you for your attention

Any Questions?

M2A has been made possible by the EU's Convergence European Social Fund through the Welsh Government



Swansea University  
Prifysgol Abertawe

**EPSRC**

Engineering and Physical Sciences  
Research Council

**COATED<sup>2</sup>**  
EPSRC Centre for Doctoral Training in Industrial Functional Coatings



UNIB EYRHWYBODAID  
EUROPEAN UNION



Llywodraeth Cymru  
Welsh Government

Cronfa Gymdeithasol Ewrop  
European Social Fund